



Ray Purt Dental

NEW PRODUCT LAUNCH



FUNCTIONS

- Reduces force through the TMJ, preventing TMJ dysfunction.
- Prevents TMJ impingement and overextension.
- Comfortably supports the jaw during long dental procedures.
- Minimises post treatment jaw pain.
- Maintains a patient airway and reduces the need for treatment breaks.
- Reduces ergonomic demands of the dental assistant and dentist, decreasing the chance of work place injuries.
- An adjustable strap at the back allows for a one size fits all, meaning even children can use this device.
- Leaves both the dental assistants hands free to assist the dentist.

"A better Dental Experience"



Three Reasons to Use JawPro

- 1) Improves Treatment Experience!
- 2) Easy to clean
- 3) Reusable

INDICATIONS

- When treating any patients (dentist, prosthodontist, orthodontist etc)
- Oral Surgery
- Anyone with a history of jaw pain
- Anxious patients or patients with movement disorders
- Procedures with high force e.g extractions
- Long Procedures
- Sedated patients who can not maintain a patent airway.

EASY TO CLEAN AND USE

The JawPro is made out of high quality plastic, making it autoclavable! The JawPro's smooth glossy surface makes cleaning simple. It can withstand strong chemicals allowing it to be easily wiped down post treatment. To put JawPro on and off is a timeless, efficient process.

REUSSABLE

JawPro has been innovated and created with the highest quality material to ensure a long product life. It has been tested against repetitive use, extreme temperatures and under high pressure in laboratory settings to ensure you are getting a strong reliable product!

HOW DOES IT WORK?

The JawPro works by displacing the downward forces exerted by the Dentist through your upper chest instead of your TMJ. The strength of the device is provided from the ratchet design that allows for an easy locking system. The device is simple to put on and is flexible allowing a one size fits all approach.

REUSABLE

JawPro has been innovated and created with the highest quality material to ensure a long product life. It has been tested against repetitive use, extreme temperatures and under high pressure in laboratory settings to ensure you are getting a strong reliable product!

To watch a video all about the JawPro please click on the following link:

https://www.youtube.com/watch?v=JezX1Sr_p8E

For Information on Fitting the JawPro click on the link below

<https://www.youtube.com/watch?v=rfVdF1a0AdU>

Contact Ray Purt Dental for further queries

Email: orders@raypurdental.com.au

Web: www.raypurdental.com.au

Telephone: 08 9367 1216